
Sorbion[®] Sorbact[®]

CUTIMED® SORBION® SORBACT® SUPER-ABSORBENT DRESSING

FOR REMOVING EXUDATE AND PHYSICALLY BINDING BACTERA








SUPER-ABSORBENT DRESSING THAT ALSO MANAGES BACTERIA


SAFE REMOVAL OF BACTERIA

- Cutimed[®] Sorbion[®] Sorbact[®] featuring Sorbact[®] Technology safely and irreversibly binds bacteria, helping to reduce toxin release and replication
- Unlike other antimicrobial wound dressings, bacteria physically bind to Sorbact[®] Technology, so cell debris may be safely removed from the wound bed
- An excellent alternative option from other antimicrobial dressings, as there is no release of active agents and no known cases of bacterial resistance.









ABSORBS AND RETAINS EXUDATE

 The superabsorbent core absorbs and retains moderate to high levels of exudate reducing the risk of skin maceration, even under compression

Ordering Information

	Size	Ref-No.	Dressings/ Box
Cutimed® Sorbion® Sorbact®			
9	4 x 4 in.	7269808	10
	4 x 8 in.	7269809	10
	8 x 8 in.	7269810	10
	8 x 12 in.	7269811	10

Indications

Cutimed® Sorbion® Sorbact® is indicated for use in the management of shallow, heavily exuding, partial to full thickness wounds (including clean, colonized, contaminated, or infected wounds): chronic wounds (venous and arterial ulcers, diabetic ulcers and pressure ulcers), postoperative dehisced wounds and traumatic wounds.

Sorbact® is a registered trademark of ABIGO Medical AB. Please refer to the product label and/or package insert for full instructions on the safe use of these products.



