

Leukomed[®] Control Post-Op Dressing

HELPING TO PREVENT SURGICAL SITE INFECTIONS



The first choice for second skin

Totally Transparent Composition

- Quick detection of first signs of infection possible
- Reduce number of dressing changes



Ultra-thin, Conforming Film That's Like a "Second Skin"

- Readily adapts and moves with skin's flexion and edema through reliable, skin-friendly adhesive border
- Highly breathable and conformable helps prevent skin maceration and tension blisters
- Transparent, low-friction, matte surface
- · High wearing comfort for patients

Thin, Absorbent Hydropolymer Wound Pad

- Ensures optimal moisture level of wound and surrounding skin
- Keeps wound clean and ready for epithelialization







Patient after shoulder surgery.



Patient after shoulder surgery.



For the correct application of Leukomed® Control, please watch our application video

Ordering Information

Leukomed® Control - Post-Op Dressings			
Size	PAD SIZE	Order Code	Dressings/Box
2 x 2 ³ / ₄ in.	³/ ₄ x 1¹/ ₂ in.	7323000	10
2 ³ / ₄ x 4 in.	1½ x 2¾ in.	7323001	10
3½ x 6 in.	1 ³ / ₄ x 4 ¹ / ₂ in.	7323002	10
4 x 9½ in.	2 ³ / ₈ x 7 ⁷ / ₈ in.	7323003	5
4 x 13 ³ / ₄ in.	2 ³ / ₈ x 7 ⁷ / ₈ in.	7323005	5

Please refer to the product label and / or package insert for full instructions on the safe use of these products.

Leukoplast*,

Indications

Leukomed® Control is indicated for dry and low exuding acute wounds such as surgical incisions, superficial partial-thickness burns, skin donor sites, lacerations and abrasions. **Recommendation:** It is especially suited for the visual control of wounds on predisposed patients with an increased risk of infection.



